


# GlucoSure™ VIVO


## Blood Glucose Monitoring System



Autocode



Fast Result



HCT Correction  
30-55%



Alternate Site  
Testing



Voice Guidance






PC Link

The GlucoSure VIVO Blood Glucose Monitoring System provides an easy and precise way to measure blood glucose using fresh capillary whole blood taken from fingertips, palms, or forearms. The GlucoSure VIVO has a large screen and Voice Guidance\* for easy understanding of operation and testing results.

**APEXBIO™**

# GlucoSure™ VIVO Blood Glucose Monitoring System

## More Testing Sites



## Easy to use



1. Easy Insertion



2. For Fingertip Testing



3. Apply Sample



4. Results in 6 Seconds

## Specifications

Measurement Range: 20-600 mg/dL

Test Time: 6 Seconds

Sample Volume:  $\geq 0.7 \mu\text{L}$

Hematocrit Range: 30-55%

Operating Temperature: 10°C-40°C

Memory: 300 Test Results

Data Transfer: PC Link

Battery: 2 Alkaline AAA Batteries


Dimensions: 93L x 58W x 20.5H (mm)

Weight: 79g

Voice Guidance\*

\*Please inquire about language available from your local distributor.

## Precision and Accuracy performance according to ISO15197:2013



| ISO 15197 : 2013   |           |           |
|--------------------|-----------|-----------|
| Sample < 100 mg/dL |           |           |
| ±5 mg/dL           | ±10 mg/dL | ±15 mg/dL |
| 131 / 180          | 176 / 180 | 180 / 180 |
| 73%                | 98%       | 100.0%    |
| Sample ≥100 mg/dL  |           |           |
| ±5%                | ±10%      | ±15%      |
| 271 / 420          | 395 / 420 | 418 / 420 |
| 65%                | 94%       | 99.5%     |

\* ©2026 Apex Biotechnology Corp. All rights reserved. APEXBIO and all product logos are registered trademarks of Apex Biotechnology Corp.

\* Specifications may vary according to different authorized bodies.

\* Actual content of the product is subject to change without notice.

