


GlucoSure[™] VIVO^{Link}

Blood Glucose Monitoring System



Autocode



Fast Result



Alternate Site
Testing



Voice Guidance



Bluetooth



PC Link




The GlucoSure VIVO Link Blood Glucose Monitoring System provides an easy and precise way to measure blood glucose using fresh capillary whole blood taken from fingertips, palms, or forearms. The GlucoSure VIVO Link has a large screen and Voice Guidance* for easy understanding of operation and testing results.

APEXBIO[™]

GlucoSure™ VIVO^{Link} Blood Glucose Monitoring System

The GlucoSure VIVO Link's Bluetooth wireless feature enables testing results to seamlessly transmit into a smart device and share the data with healthcare providers instantly for better understanding and management of patients' condition.

More Testing Sites



Specifications

Measurement Range: 20-600 mg/dL (1.1-33.3 mmol/L)

Coding Method: Automatic

Testing Time: 6 seconds

Sample Volume: $\geq 0.7 \mu\text{L}$

Operating Temperature: 10°C-40°C

Memory: 300 blood Test results with date and time


Data Management: 7, 14, and 30 days averaging

Additional Features: PC link, Bluetooth,
Alternative Site Testing

Voice Guidance*

*Please inquire about language available from your
local distributor.

Precision and Accuracy performance according to ISO15197:2013



Accuracy for glucose level < 100 mg/dL

Within ± 5 mg/dL	Within ± 10 mg/dL	Within ± 15 mg/dL
121/174 (69.5%)	170/174 (97.7%)	174/174 (100.0%)

Accuracy for glucose level ≥ 100 mg/dL

Within $\pm 5\%$	Within $\pm 10\%$	Within $\pm 15\%$
262/426 (61.5%)	401/426 (94.1%)	426/426 (100.0%)

Total within ± 15 mg/dL & $\pm 15\%$

600/600 (100.0%)

* ©2026 Apex Biotechnology Corp. All rights reserved. APEXBIO and all product logos are registered trademarks of Apex Biotechnology Corp.
* Specifications may vary according to different authorized bodies.
* Actual content of the product is subject to change without notice.

